

SURGICAL ESSENTIALS

Instrument Pouch

- Keeps frequently used instruments out of the way yet in surgeon's reach as needed
- Comes with three pockets for multiple instruments
- Extra long sleeve securely holds complete laparoscopic instrument without dangling or danger of dropping
- Two surfaces of aggressive double-sided tape for secure placement
- Not made with natural rubber latex


Product No.	Description	Quantity
25-303	3-Pocket, 18.5" x 9.5"	30/Cs
25-302	2-Pocket, 11.5" x 7"	30/Cs



Cautery Tip Cleaner

- A convenient debris removal system
- Allows quick and easy cleaning of instruments
- Coated scratch pad eliminates disruptive scraping of electrodes
- Simple design applicable to mono or bipolar units
- Radiopaque element ensures x-ray detectability
- Not made with natural rubber latex

Product No.	Description	Quantity
32-200	2" x 2"	100 Pks/Cs
30-200	2" x 2"	240 Pks/Cs



Instrument Wipe

- For scratch-free cleaning of ophthalmic and other delicate instruments
- Non-linting
- Soft, sponge-like material (polyvinyl alcohol)
- Highly absorbent
- Not made with natural rubber latex

Product No.	Description	Quantity
25-201	Instrument Wipe, Polyvinyl Alcohol (2 1/2" x 2 1/2")	120/Cs

